



I
N
T
R
O
D
U
C
I
N
G



Rebound
Medical & Safety

SARS-CoV-2 Rapid
Antigen Test Kits
Non-Invasive Anterior
Nasal or Saliva Swab
(German Government
approved for Self-Test!)


NEXT GENERATION

**2 in 1 Test Without any discomfort!
(Great for kids and sensitive Adults)**



OUR STOCK – READY TO SHIP

**RAPID
RESULTS**



If you or the kids dread the thought of another Lateral Flow COVID test, we have the answer! A pain and discomfort free alternative at a very sensible price, with outstanding performance and 15-minute results. (half the time of NHS deep swab tests)

Most Rapid Tests require swabs to be taken deep into the nose or throat. This next generation kit needs only a simple mouth or end of nose swab yet is one of the most sensitive and accurate kits on the market today.

Performance is Key

In a European independent test, 700 people were selected using laboratory PCR tests. 200 were confirmed positive with onset of the infection within 7 days. 500 were confirmed as COVID free. This kit accurately identified 192 of 200 confirmed positive cases (96% accuracy). This is one of the most accurate tests on the market and is far easier and more comfortable to use than most.



Cross Reactivity

You want to be sure that any positive result is not due to another virus. Tested with the top 30, high prevalence respiratory pathogens, such as influenza and Rhinovirus this kit registered no false positives.


Self-Testing vs. Professional Use

In a controlled test, 500 individuals were swabbed twice – Once by professional laboratory staff and the second time by the patients themselves. The professional vs. Layperson results varied by only 1.5% for sensitivity and were a 100% match for specificity.





**Rebound
Medical & Safety**



Contact Us




ppe@reboundeu.com



www.reboundeu.com/ppe



+44 (0) 1635 555966



99.6%

of all 500 subjects were able to understand the instructions and complete the test with no problems.

Easy to Use & Accurate

To ensure the best possible test accuracy, it is important that instructions are simple to follow. A QR linked video is also available to take you through the test process step by step.



Important Information

This WANTAI Test kit is approved in Germany for Self-Testing under BfArM authorisation. The test is also CE approved for professional use in all other European countries. Outside Germany, use either a healthcare professional or an individual that has been trained to an appropriate standard in testing techniques to use this product. Check out the QR Code above to see just how easy it is!

Ordering:

There are 2 order code options:



WJ-2901

Components	WJ-2901
Test Cassette	x1
Extraction Vial	0.3ml x1
Disposable Sterile Swab	x1
Plastic Sealable Bag	x1
Instructions for Use	x1
Operation Guide	x1

Please call for our latest Stock Prices and MOQ levels

WJ-2910

Components	WJ-2910
Test Cassette	x10
Extraction Vial	0.3ml x10
Disposable Sterile Swab	x10
Plastic Sealable Bag	x10
Instructions for Use	x1
Operation Guide	x1

Please call us for our latest Stock Prices and MOQ levels